

MICROAIRE COMPATIBLE CANNULAS



Alira's Single Port Cannula - MicroAire Compatible

It is used for extraction of fat during liposuction, surface irregularities are prevented by keeping port away from skin.



SKU	Diameter	Length	Port
ASPC-010	3mm	30cm	8mm
ASPC-011	3mm	22cm	8mm
ASPC-012	3mm	22cm	6mm
ASPC-013	2.4mm	22cm	7.5mm
ASPC-014	2.4mm	22cm	5mm
ASPC-015	2.4mm	15cm	7.5mm

Alira's Mercedes Cannula - MicroAire Compatible

It is aggressive extracting cannula for fat extraction during liposuction. Microaire compatible it provides complete 360 degrees extraction.



SKU	Diameter	Length	Port
AMC-010	5mm	30cm	10mm
AMC-011	5mm	22cm	10mm
AMC-012	4mm	30cm	10mm
AMC-013	4mm	22cm	10mm
AMC-014	3mm	30cm	8mm
AMC-015	3mm	22cm	8mm
AMC-016	3mm	15cm	8mm
AMC-017	2.4mm	15cm	7.5mm
AMC-018	4mm	40cm	10mm

Alira's Double Mercedes Cannula - MicroAire Compatible

Alira's Double Mercedes Cannula – MicroAire Compatible is used for aggressive extraction of fibrous tissues. It ensures large volume extraction.



SKU	Diameter	Length	Port
ADMC-010	5mm	30cm	10mm
ADMC-011	5mm	22cm	10mm
ADMC-012	4mm	30cm	10mm
ADMC-013	4mm	22cm	10mm

Alira's Mirrored Tri-Port II Cannula - MicroAire Compatible

Aggressive extracting cannula for providing 270 degree of extraction during liposuction procedure.



SKU	Diameter	Length	Port
AMTC-010	5mm	30cm	10mm
AMTC-011	5mm	22cm	10mm
AMTC-012	4mm	30cm	10mm
AMTC-013	4mm	22cm	10mm

Alira's Tri-Port II Cannula - MicroAire Compatible

Alira's Tri-Port II Cannula – MicroAire Compatible is extracting cannula for 270 degree extraction of fat during liposuction procedure.



SKU	Diameter	Length	Port
ATPC-010	4mm	30cm	10mm
ATPC-011	4mm	22cm	10mm
ATPC-012	3mm	30cm	8mm
ATPC-013	3mm	22cm	8mm

Alira's Tri-port III Cannula - MicroAire Compatible

It is moderately aggressive provides 90 degree extraction, Alira's Triport III Cannula – MicroAire Compatible is used for extraction of fat during liposuction.



SKU	Diameter	Length	Port
ATCM-010	4mm	30cm	10mm
ATCM-011	4mm	22cm	10mm
ATCM-012	3mm	30cm	8mm
ATCM-013	3mm	22cm	8mm

Alira's Multi-Hole Cannula - MicroAire Compatible

It is a non aggressive harvesting cannula, provides small lobules for small volume fat grafting. It assists in fat removal after laser and ultra sonic procedures.



SKU	Diameter	Length	Port
AMHC-010	3mm	30cm	1.9mm
AMHC-011	3mm	22cm	1.9mm
AMHC-012	3mm	15cm	1.9mm

Alira's Del Vecchio Track-12 Cannula - MicroAire Compatible

Alira's Del Vecchio Track-12 Cannula – MicroAire Compatible is aggressive harvesting cannula, designed for harvesting of small lobules of fat.



SKU	Diameter	Length	Port
ADVC-010	4mm	30cm	6mm
ADVC-011	3mm	22cm	6mm
ADVC-012	4mm	40cm	6mm

Alira's Helixed Tri-Port III Cannula - MicroAire Compatible

Alira's Helixed Tri-Port III Cannula – MicroAire Compatible is designed for high fibrous areas, aggressive for harvesting of large amount of fat during liposuction.



SKU	Diameter	Length	Port
AHTC-010	5mm	30cm	10mm
AHTC-011	5mm	22cm	10mm
AHTC-012	4mm	30cm	10mm
AHTC-013	4mm	22cm	10mm

Alira's Three Hole Spatula Cannula - MicroAire Compatible

It is non-aggressive cannula with blunt tip for safe extraction from sensitive areas like neck and back.

SKU	Diameter	Length	Port
ATSC-A010	3mm	15cm	1.5mm



Alira's Two Hole Spatula Cannula - MicroAire Compatible

It is non-aggressive cannula with blunt tip for safe extraction from sensitive areas like neck and back.

SKU	Diameter	Length	Port
ATSC-B011	3mm	15cm	1.5mm



Alira's Single Hole Spatula Cannula - MicroAire Compatible

Alira's Single Hole Spatula Cannula – MicroAire Compatible is short blunt tip cannula for fast harvesting of fat for neck region.

SKU	Diameter	Length	Port
ATSC-C012	2.4mm	15cm	3mm



Alira's Flared Mercedes Cannula - MicroAire Compatible

Alira's Flared Mercedes Cannula – MicroAire Compatible is ideal for debulking, dissection and pre-tunneling of fat before extraction.

SKU	Diameter	Length	Port
AFMC-010	5mm	30cm	10mm
AFMC-011	5mm	22cm	10mm
AFMC-012	4mm	30cm	9mm
AFMC-013	4mm	22cm	9mm
AFMC-014	3mm	30cm	8mm
AFMC-015	3mm	22cm	8mm



Alira's Bent Flared Mercedes Cannula - MicroAire Compatible

Alira's Bent Flared Mercedes Cannula – MicroAire Compatible breaks up the fat before extraction providing better harvesting and 360 degree extraction.

SKU	Diameter	Length	Port
ABFC-010	5mm	30cm	10mm
ABFC-011	4mm	30cm	9mm



Alira's Bent Mercedes Cannula - MicroAire Compatible

Alira's Bent Flared Mercedes Cannula – MicroAire Compatible breaks up the fat before extraction providing better harvesting and 360 degree extraction.

SKU	Diameter	Length	Port
ABMC-010	4mm	30cm	10mm
ABMC-011	3mm	30cm	8mm



Alira's Turbo Double Mercedes Cannula - MicroAire Compatible

Alira's Turbo Double Mercedes Cannula – MicroAire Compatible allows more air to be vent through tube during liposuction. This assists in better flow rate of the fat.

SKU	Diameter	Length	Port
ATDC-010	5mm	30cm	10mm
ATDC-011	4mm	30cm	10mm



Alira's Turbo Helixed Tri-Port III Cannula - MicroAire Compatible

Alira's Turbo Helixed Tri-Port III Cannula – MicroAire Compatible allows more air to be vent through tube during liposuction. This assists in better flow rate of the fat.

SKU	Diameter	Length	Port
ATHC-010	5mm	30cm	10mm
ATHC-011	4mm	30cm	10mm



Alira's Turbo Mercedes Cannula - MicroAire Compatible

Alira's Turbo Mercedes Cannula – MicroAire Compatible allows more air to be vent through tube during liposuction. This assists in better flow rate of the fat.

SKU	Diameter	Length	Port
ATMC-010	4mm	22cm	10mm



Alira's Straight HD Cannula - MicroAire Compatible

This microaire compatible cannula is used for extraction of superficial layer of fat during liposuction procedure.

SKU	Diameter	Length	Port
ALBC-010	4mm	30cm	10mm
ALBC-011	5mm	30cm	10mm



Alira's Left Bent HD Cannula - MicroAire Compatible

This microaire compatible cannula is used for extraction of superficial layer of fat during liposuction procedure.

SKU	Diameter	Length	Port
ASHC-010	4mm	30cm	10mm
ASHC-011	5mm	30cm	10mm



Alira's Right Bent HD Cannula - MicroAire Compatible

This microaire compatible cannula is used for extraction of superficial layer of fat during liposuction procedure.

SKU	Diameter	Length	Port
ARBC-010	4mm	30cm	10mm
ARBC-011	5mm	30cm	10mm



For More Information & Details

CONTACT US

Alira Medical Devices LLC

📞 +1 (818) 614 4592

✉️ info@alirameddevices.com

🌐 www.alirameddevices.com

📍 14622 Ventura Blvd Ste 102, Sherman Oaks, CA 91403

