

LIPOSUCTION CANNULA SET



Redefining the Art of Fat Grafting

Liposuction Cannula Set

Alira's Liposuction Cannula Set is crafted with the guidance of leading experts, ensuring optimal safety and precision for various body sculpting procedures. This comprehensive set of reusable and cell-friendly cannulas is specially designed to facilitate fine liposuction, allowing for effective contouring and reshaping of the body. Alira's Liposuction Cannula Set provides the confidence and versatility needed for all liposuction procedures. It offers the ideal cannulas for safely removing excess fat and achieving a harmonious and sculpted result, helping you meet your patients' aesthetic goals with precision and assurance.

1 Mercedes Liposuction Cannula - Three Lateral Ports

Alira's Mercedes Liposuction Cannula is a unique tool for liposuction procedures or fat harvesting. The three lateral ports aid to expedite the suction process in a variety of orientations. This circumferential three port configuration cell friendly liposuction cannula is effective for multi-plane soft tissue harvesting. In addition, this cannula is available with different sizes and working lengths depending on surgical preference. You can also choose between Luer -Lock Type connector, Super Luer Lock Type Connector, Threaded Type connector and with Fixed Power Handle as per your surgical needs.



Length: 25cm | **Diameter:** 3mm | **Qty:** 01
Length: 30cm | **Diameter:** 4mm | **Qty:** 01

Connector Type Options: Luer Lock | Ultra Luer Lock | Nouvag Threaded Hub | 60cc Toomey Hub | Power Cannula | Fixed Handle | PAL Handle

2 Basket Shaped Cannula

The three-port Basket is used in areas with skin irregularities as an equalizer, injector, and suction cannula. In addition, the cannula is available with different sizes and working lengths depending on surgical preference. You can also choose between Luer -Lock Type connector, Super Luer Lock Type Connector, Threaded Type connector and with Fixed Power Handle as per your surgical needs.



Length: 32cm | **Diameter:** 4mm | **Qty:** 01

Connector Type Options: Luer Lock | Ultra Luer Lock | Nouvag Threaded Hub | 60cc Toomey Hub | Power Cannula | Fixed Handle | PAL Handle

3 Reverse Pyramid Cannula

Alira's Reverse Pyramid Cannula is highly effective in sucking up fat tissue because it has two large upper lateral holes and one smaller center hole. For liposuction treatments or fat harvesting, this is a special tool. Additionally, this cannula comes in a variety of diameters and working lengths to suit surgical preferences.



Length: 25cm | **Diameter:** 3mm | **Qty:** 01

Length: 30cm | **Diameter:** 4mm | **Qty:** 01

Connector Type Options: Luer Lock | Ultra Luer Lock | Nouvag Threaded Hub | 60cc Toomey Hub | Power Cannula | Fixed Handle | PAL Handle

4 Spatula Cannula

An exceptional spatula cannula with a single port for liposuction, a fixed, lightweight handle, and a shaft made of German stainless steel is reusable.



Length: 15cm | **Diameter:** 3mm | **Qty:** 01

Connector Type Options: Luer Lock | Ultra Luer Lock | Nouvag Threaded Hub | 60cc Toomey Hub | Power Cannula | Fixed Handle | PAL Handle




WHY CHOOSE ALIRA?

Enhanced Precision: You can achieve precise fat graft placement, resulting in consistent and desirable cosmetic outcomes.

Minimized Patient Trauma: Alira's instruments minimize tissue trauma, leading to faster healing and reduced post-operative discomfort for patients.

Reduced Surgeon Fatigue: Ergonomic handles improve your comfort during procedures, reducing the risk of fatigue and enhancing surgical control.

Biocompatible Materials: All materials used are biocompatible, ensuring compatibility with the human body and reducing the risk of adverse reactions.



CORE AREAS OF EXPERTISE

Manufacturing
Extensive Range of
High Quality
Surgical Instruments

Instruments
Customization

Instruments
Repairing Services

Contract
Manufacturing

CONTACT US

ALIRA MEDICAL DEVICES LLC

+1 (818) 614 4592 info@alirameditives.com www.alirameditives.com

14622 Ventura Blvd Ste 102, #712 Sherman Oaks, CA 91403

