


# FHD 180B Endoscopic Camera System



## **DISCLAIMER & CONFIDENTIALITY**

This catalog is confidential and for your information only and is not intended as an offer, or a solicitation of an offer, to buy or sell any investment or other specific product. Although all information and opinions expressed in this document were obtained from sources believed to be reliable and in good faith, no warranty, express or implied, is made as to its accuracy or completeness. All information and opinions are subject to change without prior notice. Any information herein are not binding upon any party and cannot constitute any evidence in any court of law.


This catalog may not be reproduced in any form or copies circulated without prior written authority of ALIRA MEDICAL DEVICES LLC

## About us

Alira Medical Devices LLC is a California based Medical Devices' provider that offers a wide range of world-class medical devices and surgical instruments in the fields of Plastic & Reconstructive Surgery, Breast Augmentation, Liposuction, Gynecology, Laparoscopy, Urology, ENT, Ophthalmology and Dentistry to facilitate surgical procedures in a sophisticated manner.

At Alira, we focus on finer details and are aware of every part of the process, and especially, have a willingness to create value for the customer.

We have created an inspiring culture that empowers our people and provides a comfortable, functional, and fun environment that result in a team that works smarter, collaborates, and enjoys what they do.



## Quality Focus

Alira Medical Devices is committed to provide top Quality medical devices and surgical instruments to its customers by using a robust Quality Management System (QMS) derived from basic principles laid down by international guidelines e.g. GMP, ISO 13485:2016 & ISO 9001:2016 while meeting statutory/regulatory requirements.

At Alira Medical Devices, Quality begins right from the input materials. All vendors and materials are evaluated for quality & compliance attributes and regulatory requirements through a comprehensive vendor qualification and management program.

To ensure the quality of products, Alira Medical Devices is committed to develop all products following the QBD (Quality by Design) principles and a risk-based scientific approach or as per the regulatory requirements to achieve the quality targets.

Alira Medical Devices is committed to continuously upgrading the equipment and premises to fulfill the compliance requirements and deploy qualified, validated technologies and systems.

Alira Medical Devices promotes the culture of "Speak-up for Quality", shared learning, transparency, and data integrity. Also ensuring a timely and effective communication and escalation process exists to raise quality issues to the appropriate levels of management, advocating continual improvement

FHD Endoscopic Camera System



ENT
Ophthalmology
Gynecology
Laparoscopy
Hysteroscopy
Urology
Arthroscopy
General surgery

FHD Camera Control Unit



Output Adapter	HDMIx1 DVIx1 AVx1
Storage Adapter	USB2.0x1
Video Resolution	1920x1080
Power input	AC100V-240V,50/60Hz
Weight[kg]	3kg

FHD Camera Head



Adapter	The image resolution of adapter reaches 250LP/mm; The center resolution and edge resolution of 4K zoom interface reach 250LP/mm
Cable Size	\$5mmx3m
Weight[g]	176g

LED Cold light source




Output Adapter	Universal Adapter of light guide cables
Brightness	1500lm
Power input	AC 100V-240V,50/60Hz,400VA
Weight [kg]	2.5kg

Light guide cables



Length	3000mm
Effective Core Diameter	4.8/3mm
Color	Grey

4K Endoscopic Medical Monitor



Display Size	27"
ResolutionAspect Ratio	3840X2160
Input	HDMI 2.0 X2
Power supply	AC100V-240V 50/60HZ

# Quality Instruments Enhanced Practices Better **Lives!**



**Tell:** (818) 614 4592

**Email:** [info@alirameddevices.com](mailto:info@alirameddevices.com)

**Web** [www.alirameddevices.com](http://www.alirameddevices.com)



**Address:** 14622 Ventura Blvd Ste 102,  
#712 Sherman Oaks, CA 91403